






# BD FACSDiva 8.0.1





Tube: HeLaSAM

Population	#Events	%Parent	%Total
All Events	50,000	####	100.0
scatter	36,453	72.9	72.9
S1	29,133	79.9	58.3
s2	27,143	93.2	54.3
P1	1,582	5.8	3.2
P2	25,300	93.2	50.6
High	279	1.1	0.6
Low	268	1.1	0.5
P3	186	0.7	0.4
P4	232	0.9	0.5

BD FACSDiva 8.0.1





**Reset Threshold for your sample.**  
**Check your Window Extension is set 2.0**  
**Lower scatter voltage for your target population.**  
**Check the presion mode in your sort layout**